



Advanced thinking

BiPAP autoSV Advanced further simplifies management of complex sleep apnea.

Applying advanced thinking can lead to simpler solutions for your patients ... and for you.

BiPAP autoSV Advanced blends proven technologies to make treatment and management of complicated sleep apnea cases easier and more effective. The first thing you'll notice is how auto EPAP plus an advanced back-up rate make titration simpler. Then how our enhanced algorithms and sophisticated servo ventilation continually monitor and adjust to changing patient therapy needs. To round it out, our renowned Bi-Flex technology assures patient comfort and promotes long-term compliance. It all adds up to advanced thinking that helps you make the complex quite a bit simpler.

PHILIPS

RESPIRONICS

Product specifications	
Pressure range	IPAP: 4 cm H ₂ O to 30 cm H ₂ O EPAP: 4 cm H ₂ O to 25 cm H ₂ O
Breath rate	Auto, 4–30 bpm, off
Dimensions	24.8 cm (l) x 16.8 cm (w) x 11.2 cm (h)
Weight	1.8 Kg
Filters	Reusable pollen and optional ultra-fine
Device set-up	On-board keypad and Encore SmartCard
Data storage capacity	6 months @ 3 usage sessions/day; 7 days in-depth data analysis
Data storage content	Date, time, duration at pressure, short-term detailed events, long-term trends
Altitude compensation	Automatic
Electrical requirements	100 – 240 VAC 12 VDC (when operated with the external DC power adapter accessory)
DC power	Direct-connect cord
Warranty	2 years
Miscellaneous	Integrated heated humidification, lighted keys on pad, automatic on/off, mask off alert
Ordering information	
BiPAP autoSV Advanced systems	Part number
BiPAP autoSV Advanced with Encore SmartCard (device only)	1042900
BiPAP autoSV Advanced with Encore SmartCard (core package; includes device and heated humidifier)	1042901
Accessories	Part number
Pollen filters (reusable; 2 per pack)	1005964
Ultra-fine filters (disposable; 2 per pack)	1005945
Filter cap	1008479
Carrying case	1005965
Reusable patient tubing (1.8 m, gray)	622038
DC power adapter	1012975
Encore software and SmartCard	Part number
Encore Pro software	1054785
Encore USB SmartCard reader/writer	1014352
Pack of Encore SmartCards (10 cards)	1003304
Humidifiers	Part number
Heated humidifier	1005792
Heated humidifier chamber (replacement)	1035162

Respironics, autoSV, Bi-Flex, BiPAP, Encore, and Encore SmartCard are trademarks of Respironics, Inc. and its affiliates.



©2009 Koninklijke Philips Electronics N.V. All rights are reserved.
Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.
CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

Hoehc KW X/XX/XX MCI 4102401 PN 1063371

Philips Healthcare is part of Royal Philips Electronics

Europe, Middle East, Africa:

+33 1 47 52 30 00

Asia Pacific:

+65 6298 1088

Philips Respironics

1010 Murry Ridge Lane

Murrysville, PA 15668

+1 724 387 4000

+1 800 345 6443

www.philips.com/respironics